

ST 650 EQUIPMENT DISINFECTANT

- Disinfects areas on equipment missed by staff
- Reduce clinical time spent disinfecting equipment
- 360-degree disinfection to your entire device
- ULTRAVIOLET-C LAMPS: ST-600 features 8 UV-C lamps that eliminate up to 99.9% of bacteria in a few minutes
- SAFE AND SECURE TO USE: The sanitizer's light automatically shuts off as soon as you open the door. This ensures user's protection from potential UV exposure
- CHEMICAL-FREE: The sanitizer operates on UV light technology that is effective in cleansing your phone without any harsh chemicals
- PORTABLE: Thanks to its sturdy wheels, you can move the ST-600 to many areas



EFFECTIVE

- 4-log disinfection

FAST

- Can disinfect the surface of the objects in minutes with (99.99% disinfection Environment friendly
- Chemical-free disinfection, no residuals on the object surface.
- No collateral damage to the environment
- Ozone free (Ozone is a poisonous gas; it creates irritation and especially people with respiratory problems such as asthma should not be exposed)

SAFETY

- The product is fully compliant with the safety standard IEC60335-1: 2010
- Direct exposure to UV-C is dangerous to living beings, the chamber only starts when the door is securely closed, and the disinfection cycle is activated
- Auto power off when the chamber is open ensuring no UV-C exposure to the user

ROBUST

- Stainless-Steel chamber with sturdy wheels
- Easy to use, one-touch operation

Based on the data available from the National Emerging Infectious Diseases Laboratory (NEIDL)¹ at Boston University, which will be the subject of a forthcoming scientific publication by Boston university in a laboratory setting², Signify's UV-C light sources irradiating the surface of a material inoculated with SARS-CoV-2 (the virus that causes the COVID-19 disease) at a UV-C dose of 5mJ/cm² (exposure time 6 seconds) resulted in a 99% reduction of the SARS-CoV-2 virus present on the surface. The same study determined that a UV-C dose of 20mJ/cm² resulted in a reduction of 99.999% of SARS-CoV-2 virus on the surface (exposure time 25 seconds)

Virus	Host	Lamp Type	Fluence UV dose(mJ/cm2)						Protocol	Notes	Reference
			UV-1	UV-2	UV-3	UV-4	UV-5	UV-6			
SA-11	Monkey kidney Cell line MA104	LP	8	15	27	28			Yes		Sommer et al. 1989
	MA104 cell line	LP	20	80	140	200			Yes		Caballero et al. 2004
SA-11	MA104 cell line	LP	7	15	25				Yes		Chang et al. 1985
SA-11	MA104 cell line	LP	9	19	26	36	48		Yes		Wilson et al. 1992
SA-11	MA104 cell line	LP	7	15	23				Yes		Battigelli et al. 1993
SA-11 ATCC VR-1565 method: cell culture; assay based on CPE	MA104 cell ATCC CRL- 2378.1	LP	7	15	31 + tailing				Yes		Li et al. 2009
SA-11 ATCC VR-1565 method: RT- qPCR assay	MA104 cell ATCC CRL- 2378.1	LP	29	58	88	117 + tailing			Yes		Li et al. 2009
Human (HRV- Wa)	N/A	LP	16	24	32	40			Yes		Hu et al. 2012
SA-11	MA104 cell line	LP	10	21	32	43	52		Yes		Wilson et al. 1992
Siphoviridae	<i>E. coli</i> C	LP	1.8	3.6	5.7	7.5	9.3		Yes		Shin et al. 2005
T1											
	<i>E. coli</i> CN13	LP	N/A	N/A	N/A	13			Yes		Rodriguez et al. 2014
	<i>E. coli</i> CN13	MP	N/A	N/A	N/A	19			Yes		Rodriguez et al. 2014
T1UV											
Her 468	<i>E. coli</i> CN13 ATCC 700609	LP	N/A	8.3					Yes	Action spectrum	Beck et al. 2015
Her 468	<i>E. coli</i> CN13 ATCC 700609	Laser 254 nm	4.3	8.5	13	17			Yes	Action spectrum	Beck et al. 2015
T4											
	<i>E. coli</i> C	LP	1.1	2	3	4	6.7		Yes		Bohrerova et al. 2008
	<i>E. coli</i> C	MP	1.1	1.7	2.6	4	7		Yes		Bohrerova et al. 2008
	<i>E. coli</i> C	LP	3.6	8	13				Yes		Hu et al. 2012
ATCC 11303	N/A	LP	3.7	7.4	11	17	23	28	Yes		Timchak & Gitis 2012